

CHISON
Value Beyond Imaging



Redefine the Future of Portable Ultrasound

SonoAir

by CHISON



www.chison.com

[Ultra-portable]



2.1Kg

26mm



- The world's lightest and thinnest portable ultrasound : 2.1 kg & 26 mm
- Aluminum metal material for durable use.



Ultra-Connectable

- Maximum Probe ports in the world, Up to 4 , easy to connect



Ultra-Touchable

- 14.1 inch vibrant touch screen, Smoothly swipe, pinch, or zoom to optimize your image
- built-in trackpad allows you to achieve daily exams more effortlessly.
- Unique touch ID, Strong protection for your privacy and data security



Ultra-Flexible

- One machine, Two possibilities



- Dedicated trolley with height adjustable, small size, easy to move



Ultra-FAST

- 25s to power on.
- 3s to shut down.
- 3s to wake up.
- Latest SSD (128G)



Ultra-GREEN



- Battery life: 180mins+ for scanning
- 36hours+ for standby mode
- Higher integration ,Lower consumption

Ultra- Apporachable

- Innovative **3D** user interface,Designed for User Experience
- A more familiar way to work more efficiently



WIDE APPLICATION



◀ General Practice



◀ Anesthesiology



◀ Vascular Access



◀ Physical Therapy




◀ Critical Care



◀ Other POC






Bag



7L: 7.5MHz
Linear



3C: 3.0MHz
Convex



6MC: 6.0MHz
Micro-convex




3S: 2.5MHz
Phased Array




6E: 6.0MHz
Endocavity


PREMIUM PERFORMANCE




Cardiac, M Mode




Cardiac, B Mode




IVC, B Mode




Fetal Heart, B/BC Mode




Liver, B Mode




Gestational sac, B Mode




Shoulder, B Mode




Elbow, B Mode




Carotid, B/BC Mode



Ulnar artery, C Mode.



Ulnar artery, PW Mode



Thyroid, B Mode

CHISON

Value Beyond Imaging



— CHISON Medical Technologies Co., Ltd. —

Sales & Service Contact Address: No.3, Changjiang South Road, Xinwu District, Wuxi, Jiangsu, China 214028

TEL : 0086-510-85310937 **FAX :** 0086-510-85310726 **EMAIL :** export@chison.com.cn

We reserve the right to make changes to this catalogue without prior notice. All rights reserved.
Please contact our local dealer for the latest information.