







Shenzhen Enmind Technology Co., Ltd.

Add.: 5th Floor, Block A, Defengsheng Building, No.41 Dabao Road, Bao'an District 23, Shenzhen 518101, China
Tel.: +86 755 23316007

Fax: +86 755 22276344

Brand		ENMIND			
Image					
Function/Model		EN-V9 vet	EN-V7P	EN-V5 vet	EN-V3 vet
Screen		7" color touch screen	4.3" color touch screen	4.3" color touch screen	2.8" screen
Infusion mode	ml/h mode (include time and rate)	●	●	●	●
	Body weight mode	●	●	●	●
	Drip mode	●	●	●	●
	Dose mode	●	-	-	-
	Ramp up/down mode	○	○	-	-
	Sequence mode	-	○	○	○
	Drug library mode	-	○	○	○
	Loading dose mode	○	○	-	-
	Body surface mode	○	-	-	-
	Piggyback mode	○	-	-	-
	Cascade mode	●	○	-	-
	Intermittent Mode	○	-	-	-
	Micro mode	●	●	●	●
Flow rate range		0.1-2300ml/h (min.step 0.1ml/h)	0.01-2000ml/h (min.step 0.01ml/h)	0.1-2000ml/h (min.step 0.1ml/h)	0.1-2000ml/h (min.step 0.1ml/h)
WIFI		○	○	○	○
KVO		0.1-10ml/h	0.01-5ml/h	0.1-5ml/h	0.1-5ml/h
Occlusion alarm level		12 Level	4 Level or 12 Level (optional)	5 Level	4 Level
		150-975mmHg	225-900mmHg/150-975mmHg	50-900mmHg	150-900mmHg
Single bubble level		9 Level	7 Level	5 Level	5 Level
		10-800μL	50-800μL	50-800μL	50-800μL
Drug library (maximum)		5000	30/2000 (optional)	2000	30
History records(maximum)		5000	5000	5000	2000
Brand library (maximum)		200	30	200	200
Waterproof level		IP44	IP24	IP44	IP44
Appearance size parameter (Unit: mm)		218*164*82	235*95*120	138*131.5*90	138*131.5*90
Running time (minimum)		4hours @ 25ml/h	9hours @ 25ml/h	10hours @ 25ml/h	5hours @ 25ml/h
Battery capacity		2600mAh	2600mAh	5200mAh	2600mAh
Weight (bare machine does not contain fastening clips, handles and other accessories)		2kg	1.6kg	1.45kg	1.45kg
Compatible with cage fixation		-	-	●	●
Compliance with EN1789		-	●	●	●

Remarks: ● means standard, ○ means optional, --means no such function;







Shenzhen Enmind Technology Co., Ltd.

Add.: 5th Floor, Block A, Defengsheng Building, No.41 Dabao Road, Bao'an District 23, Shenzhen 518101, China

Tel.: +86 755 23316007

Fax: +86 755 22276344

Brand		ENMIND			
Image					
Function/Model		EN-S9 vet	EN-S7P	EN-S5D Vet	EN-S3 Vet
Screen		7" color touch screen	4.3" color touch screen	4.3" color touch screen	2.8" screen
Infusion mode	ml/h mode	●	●	●	●
	Bodyweight mode	●	●	●	●
	Dose mode	●	-	-	-
	Sequence mode	○	○	○	-
	Drug library mode	-	○	○	-
	Ramp up/down mode	○	○	○	-
	Loading-dose mode	●	○	○	●
	TIVA mode	○	○	○	-
	Body surface mode	○	-	-	-
	Intermittent mode	○	-	-	-
	Cascade mode	●	○	●	-
	Micro mode	●	●	○	-
Flow rate range		0.01-2300ml/h (min.step 0.01ml/h)	0.01-2000ml/h (min.step 0.01ml/h)	0.1-2000ml/h (min.step 0.1ml/h)	0.1-1500ml/h (min.step 0.1ml/h)
Syringe size		1-50/60ml	2/3-50/60ml	2/3-50/60ml	5-50/60ml
WIFI		○	○	○	○
KVO		0.1-10ml/h	0.1-5ml/h	0.1-5ml/h	0.1-5ml/h
Occlusion alarm level		12 Level	4/12 Level(optional)	4 Level	3 Level
		150-975mmHg	225-900mmHg/150-975mmHg	150-900mmHg, 50mmHg(optional)	150-900mmHg
Drug library (maximum)		5000	30/2000(optional)	2000	20
History records(maximum)		5000	5000	5000	5000
Brand library (maximum)		200	200	200	200
Waterproof level		IP44	IP24	IP44	IP24
Appearance size parameter (Unit: mm)		218*165*82	235*90*123	138*131.5*90	230*90*123
Running time (minimum)		5hours @ 5ml/h	12hours @ 5ml/h	5/10hours @ 5ml/h	12hours @ 5ml/h
Battery capacity		2600mAh	2600mAh	2600mAh	2600mAh
Weight (bare machine does not contain fastening clips, handles and other accessories)		2kg	1.8kg	2.7kg	1.8kg
Compliance with EN1789		-	●	●	●

Remarks: ● means standard, ○ means optional, --means no such function;